# Standard Operating Procedure Rigaku SmartLab XRD

I. Powder X-ray Diffraction (PXRD)



Yale West Campus Materials Characterization Core *ywcmatsci.yale.edu*  ESC II, Room E119C 810 West Campus Drive West Haven, CT 06516

- > Please FOLLOW the SOP strictly to keep the facility in good condition. Any explorations are strongly prohibited unless permitted by lab manager
- > **Only** use your Core USB drive on the XRD computer.
- > **NEVER** surf the web on the <u>SEM/EDS</u> computer in order to minimize the risk of the computer being hacked
- Yale West Campus MCC facility users must acknowledge MCC in their publications that rely significantly on MCC resources. The general acknowledgement for XRD should read:
   "The XRD data was taken using the Rigaku SmartLab X-ray at Yale West Campus Materials Characterization Core (MCC)."
- > The core reserves the right to use the data for core promotion

# Table of Contents

1	Introduction	. 1
2	Sample Preparation	2
3	Hardware Status Check	3
4	User Package Definition	. 4
5	Sample Loading	5
6	Optics Alignment	6
7	Sample Alignment	. 7
8	Sample Measurement	. 9
9	Checklist after Experiment	10

# **Rigaku SmartLab XRD Standard Operating Procedure**

- 1 Introduction
  - a) Instrument features:
    - > Full automated alignment under computer control
    - > A high-efficiency 2D detector (Hypix 3000) with high-count rate
    - > Cross Beam Optics (CBO) permits easy switching between focusing (BB) and parallel beam (PB) geometries without reconfiguration
    - > In-plane diffraction arm for in-plane measurements without reconfiguration
    - > High temperature (~1500 C) measurements in air, vacuum and helium
    - > Air-Tight Sample Container for contamination free measurement
  - 2) Location

Materials Characterization Core Room E119 810 West Campus Drive West Haven, CT 06516

b) Primary Staff Contact

Dr. Min Li Tel: 203-737-8270 Email: min.li@yale.edu Office: ESC II, Room E119D

Zishan Wu, Lab Assistant zishan.wu@yale.edu 203-824-5563 (cell) Office: ESC II

Yiren Zhong, Lab Assistant yiren.zhong@yale.edu 203-710-9820 (cell)

The Yale West Campus MCC Facilities are operated for the benefit of all researchers. If you encounter any problems with this facility, please contact the staff member listed above immediately. There is never a penalty for asking questions. If the equipment is not behaving exactly the way it should, contact a staff member.

**Notice**: Please **follow** strictly the **SOP** to keep the facility under good condition. We **DO NOT** recommend user explorations on program unless endorsed by core manager.

#### 2 Sample Preparation

**Note: always wear gloves** for sample preparation! Clean the sample holder with alcohols (IPO/ methanol)

1) Powder samples

Large particle samples **should be milled to below 5 microns** to minimize measurement deviations due to particle size effects.

#### Note:

- > It is highly recommended to use the **zero-background sample holder** so that the hump signal in the lower 2theta range (<20 degree) on regular glass holder can be avoided.
- > If sample needs to be placed on SS plate, make sure use small size Length limiting Incidence Slit (IS(L)). A sharp diffraction peak can be seen on SS plate.
- a) Small amount powder samples (~ 5 mg)
  - > Choose the **zero background sample holder** and fill the 5 mm x 0.2 mm well with powder sample and press the powder surface with a glass slide so the sample surface flushes with the holder surface

**Note**: to avoid electrostatic discharge (ESD) on magnetic materials, cover the powder with a **weighing paper** first.



> Alternately, to increase sample exposure area to X-ray beam as diffraction signals are more sensitive on sample exposure area than sample depth, the regular glass holder with  $20 \times 20 \times 0.2$  mm square indent (the shallowest one) can be used. Spread sample at the indent center and flat the surface with an optical slide or flat spatula.

**Note**: an amorphous hump from glass may appear if sample area is smaller than the X-ray beam size.



- > Fill the 20 mm x 20 mm square indent on the glass holder with the sample
- Flat the sample surface to flush with indent edge using an optical slide.
   Note: For multiple samples, the sample height alignment is required only for the first sample if all samples are prepared in the same way.
- 2) Film samples
  - a) The film samples can be placed either on zero background sample holder or regular glass holder. Choose right Incidence Slit sizes (IS1 and IS(L)). Normally IS(L) should be chosen as half of the sample size across the X-ray beam).
  - b) Alternately the stainless steel plate with the right height spacer can be used if the X-ray beam size is chosen smaller than the sample size.

## 3 Hardware Status Check

- Open the SmartLab Guildance software if it was closed (login: administrator, password: rigaku)
- 2) Check the Hardware Status window on the bottom left of the software:
  - a) Make sure that the "CBO selection slit" is "BB (Bragg-Brentano focusing mode)" Note: if it shows "PB (Parallel Beam mode), the Optics Alignment is required before measurement (see step 4)
  - b) Check if the X-ray is at standby (Tube voltage: 20 kV and Tube current: 2 mA).

If yes, click **Startup...** button in the **Package/Macro Measurement** window. In **Startup** window, choose "**Use everyday**" if the machine was used within 24 hours and hit "**Execute**" button.

**Note**: it will take 15 minutes for the X-ray power to resume the operation power (**10 kV** and **2 mA**).



## 4 User Package Definition

 On the Package window at the right side of software window, find Preinstalled> General Measurement > Bragg-Brentano Focusing package



2) **Double Click** on **General (Bragg-Brentano Focusing)** to open the package in the **Package/Macro Measurement** window:

😔 File Edit Tasks View Control Pro	cess										
Package/Macro Measurement	ìΧ										
General (Bragg-Brentano focusing)											
New Package Bar Save As											
	(										
General (Bragg-Brentano focusing)											
Startup  Key Execute Execute											
1 Optics Alignment (BB)											
<u>.</u>											
2 Sample Alignment (BB)											
Ţ											
3 General Measurement (BB)											

3) Click on Save As button in above window to save the package as user defined package which will appear in the User defined tab on the right side of the Guidance window. On the bottom of User defined tab click to create a user folder and move the package file into the folder by clicking.

## 5 Sample Loading

1) Press the Door Lock button on the Rigaku SmartLab instrument cabinet door and slide the door open



2) Push the glass sample holder or zero background holder into the **height reference sample plate** 

#### Note:

- Make sure the sample holder aligns to the center of the height reference plate
   Warning: large misalignment will lead to X-beam incidence on sample holder)
- > If no enough sample to cover the entire square well on the glass holder, try to extend the sample area along the longer side of holder (dashed line)
- > Close the cabinet door and hit **Door Lock** button



- 6 Optics Alignment
  - Put the "BB" slit into the CBO adapter and click on <sup>1 Optics Alignment (BB)</sup> in the Package/Macro Measurement window, hit "Execute" and follow the instruction on the popped up SmartMessage window:

Optics Alignment (BB)										
Change optics										
Current attribute :	Medium resolution parallel beam/RS									
Destination attribute :	: Bragg-Brentano focusing									
Optics alignment condition	Optics alignment conditions									
Change optics withou	t alignment									
Optics alignment name :	@BB focusing									
Print out results after alignment										
	2									
Execute	Import Export OK Cancel									

Press OK on the SmartMessage window after finish required steps.
 Note: this step may take about 5 minutes. Once finishes, the Optics Alignment (BB) window will be active again.

SmartMessage	
Hide figures	
Place the height reference sample plate on the attachment.	
height reference sample plate	Ľ
Insert Center_slit in height reference sample plate.	
Remove K-Beta filter from receiving slit box #2.	
OK	

3) Click OK on the Optics Alignment (BB) window after the optics alignment finishes

# 7 Sample Alignment

- 1) Click on <sup>2 Sample Alignment (BB)</sup> in the Package/Macro Measurement window
- 2) On the popped up **Sample alignment conditions** window:
  - a) Choose **Curved sample** (**Z scan only**) for particle samples, or choose **Flat sample** for flat samples such as thin films, single crystals.

b) Select Sample thickness (mm): 2.0 (if powder sample surface flushes with glass holder, or higher numbers (2.0 + solid sample thickness)
 Ting: for multiple samples with similar thickness or powders flushing with the

**Tips**: for multiple samples with similar thickness or powders flushing with the holder, **NO NEED** to repeat sample alignment.

c) Check the box next to "**Put a sample when the sample alignment starts**" and click **Execute** 

Sample Alignment (BB)
Sample alignment conditions
🗇 No height alignment
Ourved sample (Z scan only)
○ Flat sample
Sample thickness (mm) 2.0
Run recommended sequence      Customize conditions     Customize
☑ Put a sample when the sample alignment starts.
Put a sample every time the sample alignment starts in a repeated measurement.
Print alignment result.
?
Execute Import Export OK Cancel

d) Follow the instructions in the popped up SmartMessage window

Note: just hit OK if the step 4) Optics Alignment is skipped.



**Note**: the sample alignment takes about 2 minutes after which the **Sample Alignment** (**BB**) window will become active.

e) Hit OK on the Sample Alignment (BB) window

#### 8 Sample Measurement

- 1) Click on the General Measurement (BB) tab on the Package/Macro Measurement window and perform following steps:
  - a) Change File name and folder
  - b) Check K beta filter method
  - c) Select **1D mode**
  - d) Read Current Slits
  - e) Choose the right IS L size for samples (slit size: sample size = 1:2), typically 5.0 Warning: do not change other parameters on the same line.
  - f) Set typical Measurement conditions:
    - > Check Exec 1 box  $\rightarrow$  Continuous (Mode)  $\rightarrow$  Absolute (Range)  $\rightarrow$ 10.0000 (Start)  $\rightarrow$  90.0000 (Stop)  $\rightarrow$  0.0500 (Step)  $\rightarrow$ 10.0000 (Speed)  $\rightarrow$ 2/3 (IS)  $\rightarrow$ 20.000 (RS1)  $\rightarrow$  20.0 (RS2)  $\rightarrow \dots \rightarrow$ 40 (Voltage)  $\rightarrow$  44 (Current) Warning: the maximum power for X-ray source is 40 kV and 44 mA. Any higher numbers will cause severe damage to X-ray tube. Violation will lead to user account suspension and cost charges on user PI's account. Tips: to improve signal/noise ratio, users can choose smaller steps, slower speed, smaller IS, RS slit sizes or smaller Current.
  - g) Check "Drive the 4 axes to the current zero positions after measurement completed" on the bottom the General Measurement (BB) window.
  - h) Click on "Execute" to start measurement
  - i) Follow the instruction on SmartMessage windows if appear and hit OK

General N	leasurement (Bl	3)															x
Save m	easurement data																
File na	me : C:\XF	RD Use	ers\Zishan Wu\2	0160724\Sam	oleNo1.ras						•						
Sample	e name :																
Memo	:																
Manual	exchange slit cor	ndition	s					Monoc	hromatiz	ation				Da	ta acqui:	sition m	ode
Soll (	er/PSC deg)	IS L (mm)	)	PSA (deg)	Soller (deg)			🔘 K be	eta filter r	method	) Diffrac	ted beam mono	ochromator meth	od	1D	•	·
5.0	▼ 10.0	•	▼ Oper	n <b>-</b>	5.0	▼ Rea	d current s	lits									
Measur	ement conditions																
Exec.	Scan axis		Mode	Range	Start	Stop	Step	Speed	IS	RS1	RS2	Attenuator	Comment	Options	Voltage	Current	t
		_	_		(deg)	(deg)	(deg) i		deg 🔻	lmm 🔻	mm		_	_	(KV)	(mA)	
1 🔽	Theta/2-Theta	-	Continuous 🔻	Absolute 🔻	10.0000	90.0000	0.0200	10.0000	2/3	20.000	20.1			Set	40	44	-
2 🔲	Theta/2-Theta	-	Continuous 👻	Absolute 👻	3.0000	90.0000	0.0200	4.0000	2/3	2/3	0.3	-		Set	40	30	



2) XRD spectra will be automatically saved after scan

- 3) For multiple samples, change sample and start from section 3.6)
- 9 Checklist after Experiment
  - 1) For the last user of the day, click Shutdown... in Package/Macro Measurement window to set X-ray tube power to minimum (20 kV and 2 mA)

Shutdown	×
Execute	
Shutdown conditions	
XG set: Set	•
Voltage(kV): 20 Current(mA): 2	
Execute OK Can	cel

- 2) **Remove** sample from the stage
- 3) Clean the sample holders and glass slides with clean wipes using Methanol/IPA.
- 4) Store the samples holders in assigned organizer box.

- 5) **Use ONLY** the core USB flash drive to transfer data from XRD computer to the workstation in the core, and then use either your own USB flash drive or internet to retrieve data.
- 6) **Sign out** on the **logbook** and report any problems.
- 7) Keep the SmartLab Guidance software on!